

Certification by...

- ISO 9001 Certified
- PL Insurance
- Validation(IQ&OQ) Service



Innovative Technology of Ultra-Low Temperature Freezer

WiseCryo[®] WUF-Series



WUF Models -86 °C



Sales & Basic A/S:
By the Nearest Official Distributors.

Ultra-Low Temperature Freezer

Innovative SFS™ & (Simplified Freezing)




• Triple-layer Gasket for perfect insulation.



• Automatic Vacuum Dissolver making door opening easy



• Robust & Flexible designed Power Door Handle with Locking function.




• Vacuumed insulator and inner doors ensure defrosting around outer door

• Filter-Free Mechanism significantly lowers heat generated by the compressor. Block Condenser (without pins) realize filter free




WUF - 400

Comparison of Inside formation : SFS™ New technology VS competition's products.



► Complicate competition's products.



► Innovative SFS™

- Simplified Freezing System
- Much Simpler
- Oil separator free, Expansion tank free
- Block condenser against typical pin-type condenser - Dust free, Filter free

Wisdom WUF Models

