

RedLINE Etüv



By Binder redLINE ETUV

Derecelendirme:Henüz Derecelendirilmemi?
[Bu ürün hakkında soru sorun](#)

Açıklama

Technical specification:

Size (liter): 53

Temperature Range (°C): 7 above ambient to 220

Uniformity at 150 °C (± C): ± 3.6

Accuracy at 150 °C (± C): ± 0.4

Heating up time to 150 °C (min): 45

Recovery time to 150°C (min): 18

Nominal Voltage 50/60Hz (V): 230 V-CE certificate version / Plug: EU

Nominal Power (W): 800